LUXAMED®

LUXAMED® Diagnostikleuchte

Deutsch	2
English	7
Français	12
Italiano	17
Español	22
Polski	27
Português	32
Suomi	37
Svenska	42

DE

EN

FR

IT

ES

PL

PT

SU

SV

LUXAMED®

Hersteller: LUXAMED GmbH & Co. KG *

Daniel-Weil-Str.3 * 89143 Blaubeuren * Germany

Telefon: +49(0)7344 92905 0 * Fax: +49(0)7344 92905 10

E-Mail: info@luxamed.de * Internet: www.luxamed.de

- DEUTSCH -

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NE

- ENGLISH -

To ensure proper use of this LUXAMED® product, please read this user manual carefully and retain for future reference.

LUXAMED® Penlight

Intended Purpose:

The LUXAMED Penlight is intended for use for used as a throat illuminator or pupillary reflex diagnosis.

Use the device in accordance with its designated use and for its intended purpose only.

Warning:

- ⚠ Before each application, ensure that the product is externally undamaged.
- ⚠ Use only LUXAMED replacement parts and accessories. The use of third-party spare parts or accessories can adversely affect safety.
- ⚠ In order not to exceed the permissible radiation exposure, pupillary reflex diagnosis should not be performed for more than 60 seconds.
- ⚠ Always remove the batteries before cleaning / disinfecting. Ensure that no liquids enter the housing.
- ⚠ Do not service the diagnostic penlights during use on the patient.
- ⚠ Be aware of all the significant risks of reciprocal interference that may occur when using the diagnostic penlights for specific examinations or treatments.

 To avoid or minimize electromagnetic interference between the diagnostic penlight and other devices, please observe the instructions of the other devices.

 All serious incidents associated with the product must be reported to the manufacturer and to the competent authority of the Member State where the user and / or the patient is resident.

Application:

The lamp is put into operation by pressing the clip switch on the side of the Penlight. The lamp is switched off again by releasing the clip.

Replacing the batteries is done by unscrewing the lamp housing and inserting the batteries, negative pole first, into the sleeve (2 x AAA alkaline batteries with nominal voltage 1.5 V).

Blue filter (optional accessory)

With installing a blue filter on top of the LUXAMED® diagnostic light, only the blue spectral parts of the light are used.

Cleaning:

If necessary, the product can be cleaned using a damp cloth. (for example, slightly alkaline or neutral pH detergent).

Laryngeal mirror:

1. Rinse thoroughly. After disinfecting, rinse thoroughly with water. If there is stubborn dirt, such as phosphate cements, the use of "Solvitan" from VOCO is recommended.

2. Dry the laryngeal mirror very carefully.

Disinfection:

For disinfection and disinfecting cleaning of the outer surfaces the product can be applied wiping with disinfectant. For this purpose disinfectants may be used which are suitable for metal parts of medical devices. Spraying and dipping disinfection are not admitted.

Laryngeal mirror:

Steep a contaminated laryngeal mirror in a suitable disinfectant solution (such as "ID212 forte", "BIB forte", "Lysetol FF").

The following are crucial:

- Dilute the disinfectant in accordance with the manufacturer's instructions
- Steep the item for the time specified by the manufacturer's instructions

Sterilization:

Not approved for sterilization.

It is permitted to autoclave laryngeal mirrors. Place only dry laryngeal mirrors in the autoclave. No sterilization and storage in foil.

Disposal:


The product must be taken to a separate waste collection site for electrical and electronic equipment.

**Storage, transport and operating conditions:**

Ambient temperature between -10° and $+40^{\circ}$

Rel. humidity between 10% and 90%.

Symbol key

	Manufacturer
	Date of manufacture
	Observe the user manual
	Caution
	In conformity with EU Regulations.
	The device is a medical product.
	The product must be taken to a separate waste collection site for electrical and electronic equipment.
	Conforms with UK regulations
	Authorised agent for medical products in Switzerland: Inpera GmbH Thunstrasse 64 CH-3110 Münsingen, Switzerland

Guarantee Conditions.

This product is covered by a factory warranty for 2 years from the date of purchase. The date of purchase must be proven by the invoice. During the guarantee period, defects resulting from material or manufacturing faults will be repaired free of charge. The guarantee service shall not extend the guarantee period for the entire device, but only for the replaced components.

Excluded from the warranty is natural wear and tear due to normal use, transport damage and all damage that can be traced back to improper handling (for example, non-observance of the instructions for use) or to tampering with the device by unauthorized persons. The guarantee shall not be grounds for any claims against us for damages.

The statutory rights of the purchaser to claims for defects pursuant to § 437 BGB (German Civil Code) shall not be restricted.

In the event of a claim under the guarantee, the device together with the original proof of purchase must be sent to the manufacturer:

LUXAMED GmbH & Co. KG
Daniel-Weil-Straße 3
89143 Blaubeuren
GERMANY

Customer Service:

Guarantee and repair work must be carried out by trained and authorized personnel. The device must not be modified without the permission of the manufacturer. Send the device carefully packed and with properly paid postage to your authorized specialist dealer or directly to:

LUXAMED GmbH & Co. KG
Abteilung Service
Daniel-Weil-Straße 3
89143 Blaubeuren
GERMANY
Tel.: + 49 (0) 7344 92905-0

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Auktoriserat ombud för medicinska produkter i Schweiz: Inpera GmbH Thunstrasse 64 CH-3110 Münsingen, Schweiz
Advena Ltd