# LUXAMED®

**LUXAMED® Laryngoskope**

**LUXAMED® Laryngoscopes**

**Laryngoscopes LUXAMED®**

**Laringoscopio LUXAMED®**

**Laringoscopio LUXAMED®**

**LUXAMED® Laryngoskopy**

**Laringoscópios LUXAMED®**

Deutsch.....	2
English.....	12
Français.....	22
Italiano.....	32
Español.....	42
Polski.....	52
Portuguese.....	62
EMC Guidelines.....	72

# LUXAMED®

**Hersteller:** LUXAMED GmbH & Co. KG \*

Daniel-Weil-Str.3 \* 89143 Blaubeuren \* Germany

Telefon: +49(0)7344 92905 0 \* Fax: +49(0)7344 92905 10

E-Mail: info@luxamed.de \* Internet: www.luxamed.de

- DEUTSCH -

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- ENGLISH -

To ensure proper use of this LUXAMED® product, please read this user manual carefully and retain for future reference.

## **LUXAMED® Laryngoscopes**

### **Available models:**

E1.316.012 Halogen Small 2.5 V  
E1.317.012 Halogen Medium 2.5 V  
E1.318.012 Halogen Short 2.5 V  
E1.416.012 LED Small 2.5 V  
E1.417.012 LED Medium 2.5 V  
E1.418.012 LED Short 2.5 V  
E1.427.012 LED MEDIUM Charge 3.7 V  
E1.516.012 LED ECO Small 2.5 V  
E1.517.012 LED ECO Medium 2.5 V  
E1.518.012 LED ECO Short 2.5 V

### **Purpose:**

The LUXAMED® laryngoscope is a medical product to support intubation. It is used to view the larynx or the glottis.

**Side effects:** none

**Contraindications:** none

## **Warnings:**

⚠ Only blades and handles with the green mark according to ISO 7376 may be connected

⚠ LUXAMED® laryngoscope handles are intended exclusively for the operation with LUXAMED® blades.

⚠ Instrument processing must be carried out by trained and qualified personnel.

⚠ Repair work (excepted battery- and bulb changes) may only be carried out by qualified personnel. In the event of a defective rechargeable battery, please send the charging module to LUXAMED® or an authorized dealer with technical service.

⚠ LUXAMED® laryngoscopes may not be used in MRT and are not suitable for operation in explosive areas.

⚠ The instruments must be reprocessed by trained personnel.

⚠ Always remove the lamps, batteries and the charging module before cleaning and further processing steps. The charging module cannot be processed. Risk of explosion!

⚠ These instruments may only be treated with other high-quality devices to avoid flash rust of other devices.

⚠ If necessary, the charging module cover can be wiped, but no moisture may enter.

⚠ Please make sure, after each cleaning, that the plastic plate in the handle sleeve will be screwed in clockwise up to the stop.

⚠ Make sure that the head part is sufficiently screwed onto the handle (tight enough) to get a contact through the spatula fitting, contact travel according standard ISO 7376.

⚠ Make sure that the contact quality between the laryngoscope and the charging plug is adequate during charging.

- ⚠ Do not charge the device unattended.
- ⚠ When using the charging module, a rechargeable lithium-ion battery is included. Please make sure that the charging module
  - is not short-circuited
  - is not mechanically damaged
  - is not heated above the permissible temperature (see storage / transport conditions).
  - is kept away from small children
  - is always stored in a cool and dry place
  - in the event of obvious damage (leakage of ingredients, discoloration, dents, etc.), only adequate protective measures (respiratory protection, body protection) are used.
- ⚠ All serious incidents occurred in connection with the product must be reported to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.
- ⚠ The device cannot be used while charging as the power supply is disconnected.
- ⚠ Mains plug of battery charger is considered as disconnecting device.

### **Technical specifications:**

Nominal voltage	2.5V or 3.7V
Nominal current	LED 120 mA Xenon 350 mA
Applied part	Blade
Type	BF
Protection class	Class II
Ingress protection	IP20
Environment usage:	5° C -40° C 10-75% humidity 700-1070hPa
Power supply	2×AA Alkaline or 2xC Alkaline

batteries (2.5V laryngoscope)  
1x3.7V lithium secondary battery  
(3.7V laryngoscope)

**Storage / transport conditions**

temperature -10° C to 45° C.  
rel. humidity 10% - 90%  
atmospheric pressure 700-1070 hPa

**Operating conditions:**

temperature 5° C -40° C  
rel. humidity 10-75%  
atmospheric pressure 700-1070hPa

Expected lifetime of the product:5 years

**Application:**

Before each application, make sure the battery or rechargeable battery is providing sufficient capacity and the bulb is working well. When using the charging module make sure to charge approx. 5 minutes before each application.

LUXAMED® laryngoscope handles consist of handle sleeve, head section and handle cover. The LUXAMED® laryngoscope handles with rechargeable battery differ from the handle cover due to an integrated charging electronic, charging socket (Micro-USB) and charging status LED.

The integrated light is activated automatically by inserting a LUXAMED® blade onto the LUXAMED® laryngoscope handle and locking it into place. The light is deactivated by removing the LUXAMED® blade from the LUXAMED® laryngoscope handle in reverse order.

Charging the rechargeable batterie

Only use the supplied medical USB charger GTM46101-1005-USB with adapter to charge the LUXAMED® charging module. During charging process, the LED in the handle cover is lighting, after complete charging LED switches off.

The device cannot be used while charging as the power supply is disconnected.

Rechargeable battery change: in the event of a defective charging module please return it to the manufacturer or an authorized dealer with technical service.

Battery change: unscrew the handle cover and remove the used battery from handle sleeve. Insert new batteries with the positive pole in direction head section. Re-screw handle cover and test proper use of the handle.

Bulb change: unscrew the head section from the handle sleeve. Pull out the bulb from bulb holder. Now, carefully insert new bulb up to the stop. The bulb must remain in the holder by itself. Before inserting, clean the glass bulb. Re-screw head section and test proper use of the handle.

### **Battery- / rechargeable battery:**

- 2 Alkali cells (AA/LR6) for LUXAMED® laryngoscope handles paediatrics, 2 Alkali cells (C/LR14) for LUXAMED® laryngoscope handles medium
- The LUXAMED® laryngoscope handles can be used with LUXAMED® charging module with integrated charging electronics. Charging time approx. 4 hours. For this purpose, remove batteries and battery cover and insert charging module.

### **Cleaning:**

According to RKI recommendation\*, laryngoscopes are classified into semi critical medical products in accordance with group A (without special requirements on processing). Immediately after application, rinse residues under running water residue-free.

Before cleaning and further processing steps, always remove bulbs, batteries and charging module.

The cleaning agents used must be approved for medical

devices made from stainless steel, metals, plastics and glass fibres. All parts have to be residue-free before further treatment. For mechanical cleaning, always follow the instructions of the machine- and the cleaning agent manufacturer. Ultrasonic cleaning is not permitted as the product will be damaged irreparably.

### **Disinfection:**

Disinfection or disinfecting cleaning of the outer surfaces of the handle by wipe- or immersion disinfection.

Use cleaning agents approved for medical devices made from stainless steel, metals, plastics and glass fibres. Always follow the instructions of the cleaning agent manufacturer.

Thermal disinfection up to 93° C.

According to recommendations of the RKI\*, the disinfectants must apply to efficacy class AB (with proven antibacterial, antiviral and fungicidal effectiveness) according to the definition of the RKI list.

### **Sterilization:**


According to recommendations of the RKI\*, sterilization can be carried out if necessary. For this purpose, steam sterilization in autoclaves at a temperature of 134° with 5 minute hold time and a saturated steam is possible. However, different national instructions concerning hold time may exist, which do not cause any problem to the metal parts of the devices.

*\*("Requirements to the hygiene for the preparation of medical products" (Bundesgesundheitsblatt-Gesundheitsforsch-Gesundheitsschutz 2001 44:1115-1126))*

**Disposal:**

 After processing, the product must be brought to a separate collection of electrical and electronic equipment.

## Symbol key

	Manufacturer
	Date of manufacture
	Attention to the user manual
	Caution
	In Conformity with EU guidelines
	Typ BF applied part
	Class II equipment
	The product must be brought to a separate collection of electrical and electronic equipment.
	Product It is a medical device
	Conforms with UK regulations
	Authorised agent for medical products in Switzerland: Inspira GmbH Thunstrasse 64 CH-3110 Münsingen, Switzerland

## **Warranty Conditions**

From the date of purchase we give 2 years warranty. The purchase date has to be proven by the invoice. Within the warranty period defects are eliminated free of charge. After repairs the warranty period is not extended on the whole device, but only to the replaced components. Excluded from the warranty are parts subject to normal wear and tear, transport damages and any damage caused by improper handling (e.g. non-compliance with the instructions for use). Damages due to disassembly by unauthorized persons are also excluded from warranty. No claims for damages against us are substantiated by the warranty. In the case of justified warranty claims the device has to be sent along with the original invoice to:

LUXAMED GmbH & Co. KG  
Daniel-Weil-Straße 3  
89143 Blaubeuren  
GERMANY

### **Customer service**

Warranty and repair work must be undertaken by trained and authorised personnel.

Do not modify this equipment without authorization of the manufacturer.

Send the device, carefully packaged and with sufficient postage, to your authorised dealer or directly to

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
## Symbol key

	Manufacturer
	Date of manufacture
	Attention to the user manual
	Caution
	In Conformity with EU guidelines
	Typ BF applied part
	Class II equipment
	The product must be brought to a separate collection of electrical and electronic equipment.
	Product It is a medical device
	Conforme a las normas de la Reino Unido
	Representante suizo para productos sanitarios: Inspera GmbH Thunstrasse 64 CH-3110 Münsingen

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# EMC Guidelines

Below cables information are provided for EMC reference

Cable	Max. cable length, Shielded/unshielded		Num- ber	Cable clas- sification
Charger USB cable to LUXAMED® La- ryngoscope Handle (only for E1.427.012)	1m	unshielded	1 Set	DC Power

## Important information regarding Electro Magnetic Com- patibility (EMC)

LUXAMED® laryngoscope handles with charging module need special precautions regarding EMC and put into service according to the EMC information provided in the user manual and other documents.

Special precautions need to be observed:

The LUXAMED® laryngoscope handle is only to be used with a medical approved charger, plug and USB cable;

The use of accessories and cable other than those specified, with the exception of accessories and cables sold by LUXAMED® as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY or decreased LIFESPAN of the charging module.

 Wrong accessories and problems due to EMC can cause an interruption of the essential performance of supplying light for intubation.

 LUXAMED® laryngoscope handle with charging module should not be used adjacent to or stacked with other equipment. In case adjacent or stacked use is necessary, The LUXAMED® laryngoscopes should be observed to verify normal operation in the configuration in which it will be used.

Guidance and manufacturer' declaration - electromagnetic emissions - for LUXAMED® laryngoscope handle with charging module

### Emission

Test	Limit	Electromagnetic environment - guidance
Conducted emission	CISPR 11, Group 1, Class B	Device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Radiated emission	CISPR 11, Group 1, Class B	
Harmonic current emissions	IEC 61000-3-2, Class A	Device is suitable for direct connection to the public low-voltage power supply network.

## Immunity test levels

Test	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (IEC 61000-4-2)	Contact Discharge: $\pm 8$ kV Air Discharge: $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated RF EM filed (IEC 61000-4-3)	80-2700 MHz; 1kHz AM 80 %; 10 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> <p><math>d = 1.2\sqrt{P}</math> for 80 MHz to 800 MHz  <math>d = 2.3\sqrt{P}</math> for 800 MHz to 2,7 GHz            where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p>

Proximity fields from RF wireless communications equipment (IEC 61000-4-3)	385 MHz; Pulse Modulation: 18 Hz; 27 V/m 450 MHz, FM + 5 Hz deviation: 1 kHz sine; 28 V/m 710, 745, 780 MHz; Pulse Modulation: 217 Hz; 9 V/m 810, 870, 930 MHz; Pulse Modulation: 18 Hz; 28 V/m 1720, 1845, 1970 MHz; Pulse Modulation: 217 Hz; 28 V/m 2450 MHz; Pulse Modulation: 217 Hz; 28 V/m; 5240, 5500, 5785 MHz; Pulse Modulation: 217 Hz; 9 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance 30 cm.
--	--	---

<p>Electrical fast transients / bursts (IEC 61000-4-4)</p>	<p>Power lines: 2 kV; 100 kHz repetition frequency Signal lines: 1 kV; 100 kHz repetition frequency</p>	<p>Mains power quality should be that of a typical environment.</p>
<p>Surges (IEC 61000-4-5)</p>	<p>L-N: 1kV at 0°, 90°, 180°, 270° L-PE, N-PE: 2 kV at 0°, 90°, 180°, 270°</p>	<p>Mains power quality should be that of a typical environment.</p>

<p>Conducted disturbances induced by RF fields (IEC 61000-4-6)</p>	<p>0.15-80 MHz; 1kHz AM 80 %; 3 Vrms, 6 Vrms in ISM and amateur radio band</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b>  <math>d = 1.2\sqrt{P}</math> for 150 kHz to 80MHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p>
<p>Rated power frequency magnetic fields (IEC 61000-4-8)</p>	<p>30 A/m, 50 Hz</p>	<p>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</p>

<p>Voltage dips / Voltage interruptions (IEC 61000-4-11)</p>	<p>0 % UT for 0.5 cycle at 0°,45°,90°,135°,180°,225°,270°,315°  0 % UT for 1 cycle at 0°  70 % UT for 25/30 cycles at 0°  0 % UT for 250/300 cycles 0°</p>	<p>Mains power quality should be that of a typical environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device is powered from an uninterruptible power supply or battery.</p>
--	--	--

**[www.luxamed.de](http://www.luxamed.de)**

C2-1001.F.1